

## Sterile Lubricant Jelly

For lubrication prior insertion of catheters, endoscopes, and surgical instruments. Will not harm surgical instruments, rubber or plastic. Intended for medical purposes to lubricate body orifices and to facilitate entry of diagnostic or therapeutic devices.

- Gamma sterilized.
- Water soluble.
- Non irritating to skin, tissue and mucous membranes.
- Non greasy.
- Non staining.
- Non sticky.

Item No.	Description	Pkg.
53J505-4209-000	Tube 42g	12 pc/bx, 6 bx/cs
53J505-8209-000	Tube 82g	1 pc/bx, 48 bx/cs
53J505-1009-000	Tube 113g (4 oz)	12 pc/bx, 6 bx/cs
53J505-1009-011	Tube 113g (4 oz)	1 pc/bx, 56 bx/cs
53J805-3011-000	Sachets 3g	144 pc/bx, 12 bx/cs
53J805-5001-000	Sachets 5g	144 pc/bx, 12 bx/cs
53J805-5002-000	Sachets 5g	Bulk
53J805-6006-000	Sachets 15g	36 pc/bx, 144 bx/cs
53J505-1006-000	Syringes 6g	10 pc/bx, 120 bx/cs
53J505-1001-000	Syringes 8g	10 pc/bx, 120 bx/cs
53J505-1010-000	Syringes 10g	10 pc/bx, 120 bx/cs



## Non-Sterile Lubricant Jelly

Item No.	Description	Pkg.
54J505-4209-000	Tube 42g	12 pc/bx, 6 bx/cs
54J505-8209-000	Tube 82g	1 pc/bx, 48 bx/cs
54J505-1009-000	Tube 113g (4 oz)	12 pc/bx, 6 bx/cs
54J505-1009-011	Tube 113g (4 oz)	1 pc/bx, 56 bx/cs
54J805-3011-000	Sachets 3g	144 pc/bx, 12 bx/cs
54J805-5001-000	Sachets 5g	144 pc/bx, 12 bx/cs
54J805-5002-000	Sachets 5g	Bulk
54J805-6006-000	Sachets 15g	36 pc/bx, 144 bx/cs

\*Available in selected countries, for more information please contact [info@sn-medical.com](mailto:info@sn-medical.com).

## Cervilube™

Carbomer, Paraben & CHG Free.

Validated by Hologic for use with the ThinPrep® Pap test when used as instructed.

See [www.hologic.com/thinprep](http://www.hologic.com/thinprep) for more details.

CerviLube™ Lubricating Jelly is intended for medical purposes, which is used to lubricate a body orifice to facilitate entry of a diagnostic or therapeutic device.

Advantages:

- Tested safe and effective for use with Pap Test.
- Eases insertion of diagnostic and therapeutic devices.
- Transparent and greaseless.
- Sterile/Non-sterile options.
- FDA and CE regulation.
- Complies with California Prop 65, unlike other lubricants commonly used by OB/GYN offices.

CerviLube™ Lubricating Jelly does not contain:

- Carbomers - Which impact the results of the PapTest.
- Parabens - Linked to increased risk of breast cancer.
- Chlorhexidine gluconate (CHG) - Possible skin irritation and allergic reactions.

## Sterile

Item No.	Description	Pkg.
53J506-4209-000	Tube 42g	12 pc/bx, 6 bx/cs
53J506-8209-000	Tube 82g	1 pc/bx, 48 bx/cs
53J506-8200-000	Tube 82g (3 oz)	12 pc/bx, 6 bx/cs
53J506-1009-000	Tube 113g (4 oz)	12 pc/bx, 6 bx/cs
53J806-3012-000	Sachet 3g	144 pc/bx, 12 bx/cs
53J806-5011-000	Sachet 5g	144 pc/bx, 12 bx/cs



## Non-Sterile

Item No.	Description	Pkg.
54J506-4209-000	Tube 42g	12 pc/bx, 6 bx/cs
54J506-8209-000	Tube 82g	1 pc/bx, 48 bx/cs
54J506-8200-000	Tube 82g (3 oz)	12 pc/bx, 6 bx/cs
54J506-1009-000	Tube 113g (4 oz)	12 pc/bx, 6 bx/cs
54J806-3012-000	Sachet 3g	144 pc/bx, 12 bx/cs
54J806-5011-000	Sachet 5g	144 pc/bx, 12 bx/cs



\*Available in selected countries, for more information please contact [info@sn-medical.com](mailto:info@sn-medical.com).